

Monthly AED inspection

Routine readiness check. Follow the manufacturer's instructions for your model. Initial each box when the item passes.

BUILDING

PADS EXPIRE

AED LOCATION

BATTERY EXPIRES

MODEL / SERIAL

YEAR

Check	Jan	Feb	Mar	Apr	May	Jun	Jul	Aug	Sep	Oct	Nov	Dec
Status indicator shows ready, no alarm or chirp												
Pads sealed and not expired (connected if pre-connected)												
Spare pads, adult and child, not expired												
Battery not expired												
Cabinet accessible, sign visible, cabinet alarm works												
Rescue kit: gloves, razor, scissors, CPR mask, towel												
Unit clean, no cracks or damage												